


# Applix enteralsæt til pumpe, EasyBag



Enteralsæt med dråbekammer til Applix Smart pumpen. Passer til EasyBag og alle andre ernæringsposer på markedet. Passer bedst til stationær anvendelse af Applix Smart.

## **Materiale**

PVC med TEHTM\* som blødgører. Uden ftalater. Uden latex. Frigiver ikke blødgøringsmidler. Enteralsættet er gennemsigtigt, hvilket gør det let at inspicere indholdet.

Medicinporten er placeret i god afstand fra patienten og er tilpasset oral sprøjte.

\* TEHTM = Tri-2-ethylhexyl trimellitat

**Varenummer 77 51 731**

## Fresenius Kabi, Medical Devices

Article description **APPLIX PUMP SET EASYBAG**

Article number **7751731**

Material	Spike EasyBag	MABS
	Rotating nut:	PP
	Seal for Pouch Fitting:	PE
	Drip chamber enteral	SBC
	Applix Clamp:	MABS
	Roller clamp:	PP, PE
	Tube	PVC
	T-piece	MABS
	Female connector LL	PVC
	Step connector enteral	ABS

Packaging unit container	Paper	
	Foil	PA/PE

Storage container	1/30	Carton	Status: 02/11
-------------------	------	--------	---------------

PP: Polypropylene, PVC: Polyvinylchloride, ABS: Acrylnitril-Butadien-Styrene, SBC: Styrene-Butadiene-Copolymer, MABS: Methylmethacrylat-Acrylnitrile-Butadiene-Styrene