

# EasyBag Bolusadapter



Mellemstykke med spike til perforering af EasyBag-posen. Med tragtformet åbning til udtræk af sondeernæring med Janet sprøjte.

Anvendes ved bolusindgift af sondeernæring med sprøjte.

**Varenummer: 77 51 671**

## Specification

### EasyBag Bolusadapter

Fresenius Kabi Deutschland GmbH  
Production Unit Application Devices

D - 61346 Bad Homburg v.d.H.

REF	Name
7751671	EasyBag Bolusadapter EasyBag bolus adaptor

#### Zweckbestimmung / Intended use :


Adapter für Blasenspritzen zur Portionsentnahme von enteraler Nahrung aus dem EasyBag.  
Resealable attachment for EasyBag to allow removal of a bolus of enteral nutrition via a catheter tip syringe.



	Komponenten / components	material	material
1	Trichteradapter / funnel adapter	TPU	polyurethane
2	EasyBag spike	ABS	Acrylnitrole Butadiene Styrene

	Komponenten Verpackung / components: packaging	material	material
	Folie / deep drawing film	PA/PE	polyamide / polyethylene (30/70 $\mu$ )
	Papier / cover	paper	sterile paper (70 g/m <sup>2</sup> )
	Karton / carton	Cardboard	Cardboard

#### Picture:

#### Komponenten / components:



15.07.2008 Date	 Compiled by	 Approved by
--------------------	--	---