


Enteralsæt, standard til EasyBag



Enteralsæt med dråbekammer til gravitation. Med rulleklemme for justering af hastighed. Passer til EasyBag og alle andre ernæringsposer på markedet.

Materiale

PVC med TEHTM* som blødgører. Uden ftalater. Uden latex. Frigiver ikke blødgøringsmidler. Enteralsættet er gennemsigtigt, hvilket gør det let at inspicere indholdet.

Medicinporten er placeret i god afstand fra patienten og er tilpasset oral sprøjte.

* TEHTM = Tri-2-ethylhexyl trimellitat

Varenummer 77 51 041

Fresenius Kabi, Medical Devices

Article description	APPLIX GRAVITY SET EASYBAG		
Article number	7751041		
Material	Spike EasyBag	MABS	
	Rotating nut:	PP	
	Seal for Pouch Fitting:	PE	
	Drip chamber enteral	SBC	
	Roller clamp:	PP, PE	
	Tube	PVC	
	T-piece	MABS	
	Female connector LL	PVC	
	Step connector enteral	ABS	
Packaging unit container	Paper		
	Foil	PA/PE	
Storage container	1/30	Carton	Status: 02/11

PP: Polypropylene, PVC: Polyvinylchloride, ABS: Acrylnitril-Butadien-Styrene, SBC: Styrene-Butadiene-Copolymer, MABS: Methylmethacrylat-Acrylnitrile-Butadiene-Styrene